



GEMS™ Lite

The affordable alternative to hardware receiving systems

The simplest, easiest software solution for Event receiving and TTM Pacing Follow-up

GEMS™ Lite (Global ECG Management System) is a cost effective alternative to using expensive and outdated hardware receiving systems and can be installed on most existing desktop computers (refer to minimum PC requirements).

The focus of GEMS™ Lite is to manage your information from a patient centered point of view. The software captures transmitted ECG data and stores it in a simple to use database along with the associated patient information. This is the easiest way to receive and electronically store ECG's and patient information on the market today.


GEMS™ Lite is a single user, stand alone MS Access based application that allows for future expansion into the higher-end GEMS™ suite of Arrhythmia and Pacemaker/ICD management products.

GEMS™ Lite is a very affordable alternative to expensive and outdated hardware receiving systems, yet it allows you to easily expand into our higher-end GEMS™ suite of products including enhanced Event Monitoring and Pacemaker / ICD Follow-up.

GEMS™ Lite is designed to speed up and simplify basic information gathering and reporting for Arrhythmia monitoring and TTM Pacing patients. All patient and physician information is stored locally and is displayed in a format that is both easy to view and use, even for novice computer users.

GEMS™ Lite

- **efficient and easy to use**
- **cost effective ECG receiving solution**
- **improved quality of care**
- **new scheduling ability**
- **upgradeable to GEMS™ software**



**CardioComm
Solutions, Inc.**

Key Features

Reports

- Strip Report with attached Full Disclosure
- TTM Pacing Report
- 12 Lead ECG Report

GEMS™ Lite gives you three basic reporting options and each report can be saved, stored and viewed in either Adobe Acrobat™ (PDF) or Word (RTF) format.

ECG Recording

- Single and multi-lead ECG acquisition and display
- Automatic decoding of device IDs and event markers
- Compatible with most event recorders

ECG's are classified into two types, Arrhythmia / Arrhythmia 12-Lead and TTM Pacing. GEMS™ Lite has the ability to automatically split out ECGs into individual events (device dependent) and allow editing of each individual event.

ECG Editing

- On-screen caliper & beat caliper measurements
- Speed, gain and frequency filters
- Variable sample labels and lengths

In each ECG stored in the database, you can add caliper measurements, Beat Complex calipers, mark samples and add text comments.

Scheduler Lite

- Easy to navigate 'day-timer' appointment scheduler
- Auto-scheduling of recurring appointments
- Custom appointment time lengths
- Tracking and rescheduling of cancelled or missed appointments

Scheduler Lite™ is a tool for managing schedules for Arrhythmia and TTM follow-ups. Scheduler Lite is equivalent to a traditional appointment ledger. You get the familiar look and feel of an appointment book with the convenient features of a centralized database for managing all your patient appointments.

Customer Support

Customer Support

- Unmatched 24/7/365 customer service and support
- Desktop support via Internet (WebEx™ software) Software Services
- Remote installation and training

Software Services

- Remote installation and training
- Upgrade from GEMS™ Lite to GEMS™
- Custom reports (GEMS™ only)

Technical Specs

Available Database type = Microsoft Access only

Recommended PC requirements

- Intel® Celeron® Dual-Core
- Windows 7, XP, 2000
- 1 GB RAM
- 20 GB Hard Drive space
- Sound Card
- CD Rom Drive
- Internet Explorer 6 + Internet connection (Remote Desktop Support only)
- Printer



**CardioComm
Solutions, Inc.**