

# HeartCheck™ CardiBeat

Handheld ECG Monitor  
Powered by CardioComm Solutions, Inc.

## The HeartCheck™ CardiBeat Cardiac Event Monitor




Capture Lead I or II heart rhythm recordings.

Connect to a smart device with the GEMS™ Mobile App via Bluetooth to import ECG recordings.

Transmit ECGs to a physician, clinic or the SMART monitoring ECG reading service or generate your own ECG record in PDF format to show to your physician.

Available OTC the HeartCheck™ CardiBeat allows a user to record a medical grade 30 second single lead ECG. Coupled with the GEMS™ Mobile App a user can save up to a five minute recording from the device. Specific symptom, activity and diary entry information related to the recording can be entered and sent for a review by a physician, clinic or the SMART monitoring ECG reading service.

Capable of fitting into multiple clinic and telemedicine workflows, the HeartCheck™ CardiBeat with the GEMS™ Mobile App is an ideal tool for remote patient monitoring in today's ever expanding telemedicine market.



## Key Features

FDA and Health Canada cleared

Bluetooth V 4.0

Rechargeable battery

Lead: I and II

Recording Time: 30 - 300 seconds

Frequency: 0.5 to 40 Hz

Input Range: 0.05 - 6 mV

Sampling rate: 200 Hz

Size: 87 mm (L) \* 49 mm (w) \* 7 mm (H)

Weight: 33g



**CardioComm  
Solutions Inc.**

259 Yorkland Road, Suite 200  
North York, ON Canada

For more information please contact  
CardioComm Solutions, Inc. at

[sales@cardiocommsolutions.com](mailto:sales@cardiocommsolutions.com) or +1-416-977-9425 ex. 1